



LABEL WITH WATER-SOLUBLE ADHESIVE

Convenience and functionality.

LABEL WITH WATER-SOLUBLE ADHESIVE.

Accuracy in reporting information is a fundamental requirement in the medical, laboratory and pharmaceutical fields. At the same time, it may be interesting to catch the opportunity to recovering bottles and vials that have been already filled on which a label has been incorrectly applied.

To make this operation simpler, the Aro Group has developed labels with water-soluble adhesive. These are labels whose glue loses its power when in contact with lukewarm water. Removing the label becomes easier when bottles are dipped in water. After a few seconds, new labeling can be carried out, ensuring accuracy and cost savings

[LABEL WITH WATER-SOLUBLE ADHESIVE](#) [1]





© **Copyright Aro s.p.a.**

P.IVA 01723810121

Via per Cedrate, 292 - 21044 Cavarina con Premezzo (VA) - Italia

Tel. 0039 0331/212309 - Fax 0039 0331/219050

www.gruppoaro.it - aro@gruppoaro.it